




Polypectomy Snares Standard

Used for endoscopic excision of polyps using monopolar electrosurgical unit and retrieval of foreign bodies from Gastrointestinal Tract

Sheath PTFE



Art. No.	Product name	Loop Shape	Working Channel Ø	Working length
SPSO	Polypectomy Snare	Oval	2.8 mm	220 cm
SPSH	Polypectomy Snare	Hexagonal	2.8 mm	220 cm
SPSC	Polypectomy Snare	Crescent	2.8 mm	220 cm
SPSD	Polypectomy Snare	Duckbill	2.8 mm	220 cm
SPSM	Polypectomy Snare	Monofilament	2.8 mm	220 cm