

Polypectomy Snares "rotatable"

This oval snare can be rotated, so that the snare can be put into the correct position to excise the polyp properly

Operation through endoscope using a monopolar electrosurgical unit

Sheath PTFE

Monofilament for sharp cut

Multifilament for soft cut



Monofilament or Multifilament

Art. No.	Working Length	Sheath Ø	Working Channel Ø	Basket Ø
WS-2423DT-15R	230 cm	2.4 mm	2.8 mm	15 mm
WS-2423DT-20R	230 cm	2.4 mm	2.8 mm	20 mm
WS-2423DT-25R	230 cm	2.4 mm	2.8 mm	25 mm
WS-S2423DT-15R	230 cm	2.4 mm	2.8 mm	15 mm
WS-S2423DT-20R	230 cm	2.4 mm	2.8 mm	20 mm
WS-S2423DT-25R	230 cm	2.4 mm	2.8 mm	25 mm

Working Length 1650 mm available – Loop Diameter 35 mm available