


# Cytology Brush

**Used endoscopically to obtain mucosal tissue samples from Gastrointestinal Tract for microscopic-histopathological examination**



| Art. No. compatibility | Product Name               | Working Channel Ø | Length | Guidewire |
|------------------------|----------------------------|-------------------|--------|-----------|
| CY                     | Cytology Brush             | 2.8 mm            | 230 cm | 0.032"    |
| CY-DL                  | Wire Guided Cytology Brush | 2.8 mm            | 230 cm | 0.032"    |