


Biopsy Forceps

with PTFE sheath

Used endoscopically to obtain mucosal tissue samples from Gastrointestinal Tract for microscopic-histopathological examination

For easy insertion and sliding withdrawal without residues



Art. No.	Type	Working Length	Sheath Ø	Working Channel Ø
WS-2415BWH	with spike	150 cm	2.4 mm	2.8 mm
WS-2422BWH	with spike	220 cm	2.4 mm	2.8 mm
WS-2415BTH	without spike	150 cm	2.4 mm	2.8 mm
WS-2422BTH	without spike	220 cm	2.4 mm	2.8 mm