


NEG Titanium Guidewire

Nitinol

Used to gain access through difficult biliary and pancreatic strictures. To facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.

Tip straight or angled

Hydrophilic coated



Flexible tip

Art. No.	Product name	Length	Guidewire Ø	
NEG-032-S	NEG Guidewire	450 cm	0.032"	straight tip
NEG-032-A	NEG Guidewire	450 cm	0.032"	angled tip
NEG-035-S	NEG Guidewire	450 cm	0.035"	straight tip
NEG-035-A	NEG Guidewire	450 cm	0.035"	angled tip

J-tip available / wire Ø 0.025" and 0.038" available