


Hydra NEG Titanium Guidewire

Nitinol

Used to gain access through difficult biliary and pancreatic strictures. To facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.

Tip straight, super flexible – with **Extra Filiform** Wire

Hydrophilic coated



Super flexible tip

Art. No.	Product name	Length	Guidewire Ø
HYNEG-32	Hydra NEG Guidewire	450 cm	0.032"
HYNEG-35	Hydra NEG Guidewire	450 cm	0.035"
HYNEG-3226	Hydra NEG Guidewire	260 cm	0.032"
HYNEG-3526	Hydra NEG Guidewire	260 cm	0.035"

Tip angled and

Length 150 cm available (CKS)