

# Biliary Dilator

Used to dilate papilla or biliary strictures through the scope

Radiopaque PTFE material



Tip 10 Fr



Tip 7 Fr



Art. No.	Product name	Dilator Ø	Tip Ø	Working Channel Ø	Tapering Length	Usable Length
S-BD-7	Biliary Dilator	7 Fr	4 Fr	3.2 mm	3 cm	200 cm
S-BD-8.5	Biliary Dilator	8.5 Fr	5.5 Fr	4.2 mm	3 cm	200 cm
S-BD-10	Biliary Dilator	10 Fr	6 Fr	4.2 mm	3 cm	200 cm
S-BD-11.5	Biliary Dilator	11.5 Fr	8.5 Fr	4.2 mm	3 cm	200 cm