


# Biliary Balloon Dilator

**Used to dilate biliary strictures through the scope**

Shaft diameter: 2.5 cm  
Guidewire compatibility: 0.035"  
Atraumatic Tip



Art. No.	Product name	Balloon Ø	Balloon length	Working length	Burst Pressure
S-BBD-0640	Balloon Dilator	6 mm	40 mm	200 cm	3 ATM
S-BBD-0840	Balloon Dilator	8 mm	40 mm	200 cm	3 ATM
S-BBD-1040	Balloon Dilator	10 mm	40 mm	200 cm	3 ATM
S-BBD-1260	Balloon Dilator	12 mm	60 mm	200 cm	3 ATM
S-BBD-1460	Balloon Dilator	14 mm	60 mm	200 cm	3 ATM