


ERCP Cannula

Used endoscopically to opacify the biliary duct with radiopaque contrast material for fluoroscopic view, and to get easy access into the papilla.

Required guidewire 0.035"



Art. No.	Product name	Tip Ø	Working length	Working Channel Ø
SERC-T	Single Lumen ERCP Cannula Tapered Tip	5.5 Fr	200 cm	2.8 mm
SERC-B	Single Lumen ERCP Cannula Ball Tip	5.5 Fr	200 cm	2.8 mm
SERC-DL	Double Lumen ERCP Cannula Tapered Tip	6 Fr	200 cm	2.8 mm