


Self Expandable Titanium Biliary Stent

Used for stenosis or obstructed Biliary Tract caused by benign or malignant changes.

With Introducer System
Radiotherapy compatible

Uncovered stent

Stent introducer system diameter – 8 Fr
Compatible with 0.038" guidewire
Compatible with instrument channel size 3.2 mm



Art. No.	Product name	Stent length	Stent Ø
S-MBS-6	Self Expandable Biliary Stent	60 mm	10 mm
S-MBS-8	Self Expandable Biliary Stent	80 mm	10 mm
S-MBS-10	Self Expandable Biliary Stent	100 mm	10 mm