


# Single Action Stent Introducer System

## for 8.5 Fr, 10 Fr and 11.5 Fr Biliary Stents

**Used for endoscopically biliary stent placement to reduce time and easier stent placement**

Radiopaque marker over guiding catheter



Art. No.	Product name	Size	Required Guidewire	Length
SSIS-8.5	Biliary Stent Introducer System	Fr 8.5	0.038"	
SSIS-10	Biliary Stent Introducer System	Fr 10	0.038"	
SSIS-11.5	Biliary Stent Introducer System	Fr 11.5	0.038"	
Supplied without Guidewire and Stent				
BSI 5	Biliary Stent Pusher only	Fr 5		170 cm
BSI 7	Biliary Stent Pusher only	Fr 7		170 cm