



Sphincterotome

with Cutting Wire

Papillotome triple lumen – wire guided

Used for endoscopic cannulation of biliary ductal system and sphincterotomy

Working length 200 cm
Cutting Wire S.S. Monofilament
Guidewire compatibility 0.035"
Tapered tip 5.5 Fr



Art. No.	Product Name	Tip length	Cutting Loop	Remarks
SSPH-SN-20	Sphincterotome Short Nose	5.0 mm	20 mm	with Insulation
SSPH-SN-25	Sphincterotome Short Nose	5.0 mm	25 mm	with Insulation
SSPH-SN-30	Sphincterotome Short Nose	5.0 mm	30 mm	with Insulation
SSPH-LN-20	Sphincterotome Long Nose	8.0 mm	20 mm	with Insulation
SSPH-LN-25	Sphincterotome Long Nose	8.0 mm	25 mm	with Insulation
SSPH-LN-30	Sphincterotome Long Nose	8.0 mm	30 mm	with Insulation