

Endolithotripter

**Used for Stone Retrieval through the Scope
and for Lithotripsy to Crush Stones in Biliary Duct**

Light weight – with detachable flexible sheath and basket
One way directional movement locking system

Advantage: No need to remove endoscope for lithotripsy, can be done through the Scope under vision C.B.D. stone lithotripsy-retrieval procedure.

- Atraumatic Bullet type basket tip for ease of insertion and procedure
- Compatible with instrument – Channel 4.2 mm
- Basket Diameter 35 mm
- Basket Shape-Hexagonal
- Flush port for Contrast



Art. No.	Product name	Sheath Ø	Sheath length
SEL-SET	EndoLithotripter Complete Set with handle	3.7 mm	210 cm
SEL-SHEATH	EndoLithotripter Flex. Sheath with Basket Filament & Assembly	3.7 mm	210 cm
SEL-B	EndoLithotripter Basket only		
SEL-S	EndoLithotripter Sheath only		

Basket sizes 25, 30 and 40 mm available