


Manual Lithotripter

Used for Common Biliary Duct Stone Retrieval Lithotripsy

Light weight with unidirectional movement
 Sheath tip: Inner cutting edges to perform lithotripsy.
 Outer tip: rounded atraumatic
 Secure and firm fitting of SS Sheath with handle to facilitate retrieval-lithotripsy procedure



Art. No.	Product name	Sheath Ø	Sheath length
SML	Manual Lithotripter	14.5 Fr	90 cm
SML-SHEATH	Sheath for Manual Lithotripter	14.5 Fr	90 cm
SML-SHEATH-718	Sheath for Manual Lithotripter	7 Fr	180 cm
SML-SHEATH-118	Sheath for Manual Lithotripter	11 Fr	180 cm
SML-H	Manual Lithotripter Handle only		