

# DIALYZÁTORY - VitaPES®

Syntetická membrána 3. generace PUREMA®



MODEL	VitaPES® 150HF	VitaPES® 170HF	VitaPES® 190HF	VitaPES® 210HF	VitaPES® 130LF	VitaPES® 150LF	VitaPES® 180LF	VitaPES® 200LF
Katalog. číslo	70104215	70108217	70105219	70106221	70107113	70101115	70102118	70103120

## Ultrafiltration coefficient (ml/h/mmHg)

	60	67	73	80	9	10	13	15
--	----	----	----	----	---	----	----	----

## Clearances: Q<sub>B</sub> 300 ml/min

UREA	271	276	279	282	234	242	251	257
CREATIN	246	254	261	268	204	214	226	234
PHOSPHATE	232	238	242	249	165	176	189	199
VITAMIN B 12	173	185	193	198	80	87	99	105
INULIN	105	115	121	123	--	--	--	--

## Technické informace:

MODEL	VitaPES® 150HF	VitaPES® 170HF	VitaPES® 190HF	VitaPES® 210HF	VitaPES® 130LF	VitaPES® 150LF	VitaPES® 180LF	VitaPES® 200LF
Surface	1,5	1,7	1,9	2,1	1,3	1,5	1,8	2,0
Priming	89	99	112	123	82	89	108	117
Units per	30 / 960	30 / 960	30 / 960	30 / 960	30 / 960	30 / 960	30 / 960	30 / 960

**Membrane** PUREMA® Polyethersulfone (Membrana GmbH, Německo)

**Sterilizacian** Electron Beam