


Ureteral Stents

Longlife Superstent


Indwelling time up to six months and more

Made of Carbothane[®] Polyurethane
minimizes incrustation of the pigtails

1 set consists of stent (Open tip)
+ pusher, with (+G) or without guidewire



New pusher with x-ray marking



Art. No.	Size (Fr/CH)	Length (cm)	Guidewire
SOT-PRO-4826	4.8	26	0.032"
SOT-PRO-4828	4.8	28	0.032"
SOT-PRO-626	6	26	0.035"
SOT-PRO-628	6	28	0.035"
SOT-PRO-726	7	26	0.038"
SOT-PRO-728	7	28	0.038"

With closed tip and other sizes available

If you like the set with guidewire please
add +G (Example SOT-PRO 626+G)