


Ureteral Stents


Hydrophilic

Standard Double J stents made of aliphatic Polyurethane, (Tecoflex[®]) hydrophilic coated for easy insertion into the ureter

1 set consists of Double J stent (open tip) and pusher, with (+G) or without guidewire



New pusher with x-ray marking



Art. No.	Size (Fr/CH)	Length (cm)	Guidewire
SOT-4826-H	4.8	26	0.032"
SOT-4828-H	4.8	28	0.032"
SOT-626-H	6	26	0.035"
SOT-628-H	6	28	0.035"
SOT-726-H	7	26	0.038"
SOT-728-H	7	28	0.038"

other sizes available

If you like the set with guidewire please add +G
(Example SOT-4826-H+G)