


## Ureter Access Sheath

Provides **ureteral dilation** and a **continous working channel** for the introduction of Endoscopes and Instruments during ureteral access procedures

The sheath protects the ureter, delicate instruments and the smaller flexible scopes during repeated instrument exchanges, reducing the potential of trauma.



Art. No.	Sheath Inner Diameter Size (Fr/CH) = mm		Sheath Length
<b>UAS 9-13</b>	9	= 3	13
<b>UAS 9-20</b>	9	= 3	20
<b>UAS 9-28</b>	9	= 3	28
<b>UAS 9-30</b>	9	= 3	30
<b>UAS 9-35</b>	9	= 3	35
<b>UAS 9-45</b>	9	= 3	45
<b>UAS 12-13</b>	12	= 4	13
<b>UAS 12-20</b>	12	= 4	20
<b>UAS 12-28</b>	12	= 4	28
<b>UAS 12-30</b>	12	= 4	30
<b>UAS 12-35</b>	12	= 4	35
<b>UAS 12-45</b>	12	= 4	45
<b>UAS 14-13</b>	14	= 4.7	13
<b>UAS 14-20</b>	14	= 4.7	20
<b>UAS 14-28</b>	14	= 4.7	28
<b>UAS 14-30</b>	14	= 4.7	30
<b>UAS 14-35</b>	14	= 4.7	35
<b>UAS 14-45</b>	14	= 4.7	45