


# Percutaneous Nephrostomy Drainage Set "PIGTAIL"



# Percutaneous Nephrostomy Drainage Set "PIGTAIL"

A convenient complete set for percutaneous placement of drainage catheter – PIGTAIL design. Soft polyurethane material provides superior biocompatibility and drainage. Highly radiopaque for precise placement.

1 Set consists of:

- PIGTAIL – Drainage Catheter 35 cm, plus straightener, usable length 30 cm
- 2-part puncture needle 1,3 mm (18 GA)
- Lunderquist Guidewire 80 cm J-Tip, fixed core, PTFE coated 0.038"
- Fascial Dilators 20 cm
- Urine Bag connector

## Percutaneous Nephrostomy Drainage Set "PIGTAL"

Catheter hydrophilic coated

Art. No.	Catheter Size (Fr/CH)	Fascial Dilators (Fr/CH)
<b>PDS-P</b>	6	5 – 6 – 7
<b>PDS-P7</b>	7	6 – 7 – 8
<b>PDS-P8</b>	8	7 – 8 – 9
<b>PDS-P9</b>	9	7 – 8 – 9 – 10
<b>PDS-P10</b>	10	6 – 8 – 10 – 12
<b>PDS-P12</b>	12	6 – 8 – 10 – 12
<b>PDS-P14</b>	14	6 – 8 – 12 – 14

### PIGTAIL Catheter only

Art. No.	Size (Fr/CH)
<b>PCN-P6</b>	6
<b>PCN-P7</b>	7
<b>PCN-P8</b>	8
<b>PCN-P9</b>	9
<b>PCN-P10</b>	10
<b>PCN-P12</b>	12
<b>PCN-P14</b>	14

### Exchange Set

Art. No.	Catheter Size (Fr/CH)
<b>PDS-PE 8</b>	8
<b>PDS-PE 10</b>	10
<b>PDS-PE 12</b>	12
<b>PDS-PE 14</b>	14

Same as drainage set without puncture needle