


Percutaneous Nephrostomy Drainage Set "MALECOT"



Percutaneous Nephrostomy Drainage Set "MALECOT"

A convenient complete set for percutaneous placement of drainage catheter – MALECOT design. Soft polyurethane material provides superior biocompatibility and drainage. Highly radiopaque for precise placement.

1 Set consists of:

- MALECOT– Drainage Catheter 35 cm, plus straightener, usable length 30 cm
- 2-part puncture needle 1,3 mm (18 GA)
- Lunderquist Guidewire 80 cm J-Tip, fixed core, PTFE coated 0.038"
- Fascial Dilators 20 cm
- Urine Bag connector

Percutaneous Nephrostomy Drainage Set "MALECOT"

Catheter hydrophilic coated

Art. No.	Catheter Size (Fr/CH)	Fascial Dilators (Fr/CH)
PDS-M10	10	6 – 8 – 10 – 12
PDS-M12	12	6 – 8 – 12 – 14
PDS-M14	14	8 – 12 – 14 – 16
PDS-M16	16	8 – 12 – 14 – 16
PDS-M18 / M20 / M22 available		

MALECOT Catheter only

Art. No.	Size (Fr/CH)
PCN-M10	10
PCN-M12	12
PCN-M14	14
PCN-M16	16
PDS-M18 / M20 / M22 available	