

Amplatz Renal Dilators

Hydrophilic coated

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.

Renal Dilators, single

Art. No.	Size (Fr/CH)	Length (cm)
RD 18	18	30
RD 20	20	30
RD 22	22	30
RD 24	24	30
RD 26	26	30
RD 28	28	30
RD 30	30	30

Fascial Dilators, single

Art. No.	Size (Fr/CH)	Length (cm)
FD 8	8	30
FD 10	10	30
FD 12	12	30
FD 14	14	30
FD 16	16	30

Amplatz Sheath

Hydrophilic coated

Using during renal dilatation to provide atraumatic working tracks after removal of the dilator



Art. No.	Size (Fr/CH)	Length (cm)
AM 24	24	17
AM 26	26	17
AM 28	28	17
AM 30	30	17