

mediform

SPECIÁLNÍ
ZDRAVOTNICKÝ
MATERIÁL

Joline®

Allevo

Kyphoplasty System



Balloon in perfection

Catheter specialists for more than 10 years

Allevo

Kyphoplasty System

Allevo Kyphoplasty System - Advanced stabilization & pain relief for spinal fractures

Joline's minimally-invasive balloon catheters provide a complete system for treating vertebral fractures. Our unique balloon design allows our balloon catheters to be precisely positioned in the vertebra, thereby creating the best conditions for stabilizing the spine. Years of experience and precision engineering combine to provide numerous advantages and guarantee users the highest degree of safety and reliability.


Joline — Close to the customer.

Basic
Kit


Smart
Kit

Best
Kit

The convenience of having the most
suitable Kit for every Kyphoplasty treatment



Color-coding makes our
instruments quick & easy to identify



S5001

One Pedicle Kit

Balloon Ø 16 mm;
Balloon Length = 10mm

Kit contains:

1 x Vertebra Access Device
1x each Guide Pin, blunt, sharp
1 x Vertebra Introducer Device
1 x Vertebra Biopsy Device
1 x Bone Drill
1 x appropriate balloon catheter
1 x Vacuum Syringe
1 x Inflation Device
3 x Vertebra Filler Device
1 x Bone cement

S5002

One Pedicle Kit

Balloon Ø 16 mm;
Balloon Length = 16mm

S5003

One Pedicle Kit

Balloon Ø 16 mm;
Balloon Length = 22mm

Best
Kit



S1001

Best Kit 10 / 16

Balloon Ø 16 mm;
Balloon length = 10 mm

Kit contains:

2 x Vertebra Access Device
each 2 x Guide Pin, blunt, sharp
1 x Vertebra Introducer Device
1 x Vertebra Introducer Cannula
1 x Vertebra Biopsy Device
1 x Bone Drill
2 x appropriate Allevo Balloon Catheter
1 x Vacuum Syringe
2 x Inflation Device
6 x Vertebra Filler Device
1 x Bone cement

S1002

Best Kit 16 / 16

Balloon Ø 16 mm;
Balloon length = 16 mm

S1003

Best Kit 22 / 16

Balloon Ø 16 mm;
Balloon length = 22 mm

mediform

SPECIÁLNÍ
ZDRAVOTNICKÝ
MATERIÁL

MEDIFORM, spol. s r.o., Olšová 1, CZ Brno 637 00
TELEFON: +420 541 221 392 / FAX: +420 541 220 736
E-MAIL: objednavky@mediform.cz / www.mediform.cz